Python Documentation

version

June 30, 2020

Contents

Welcome to moseq2-extract's documentation!	1
moseq2_extract package	1
CLI Module	1
moseq2-extract	1
aggregate-results	1
convert-raw-to-avi	1
copy-slice	1
download-flip-file	2
extract	2
find-roi	4
generate-config	5
generate-index	5
version	6
GUI Module	6
General Utilities Module	8
Subpackages	13
moseq2_extract.extract package	13
Extract - Extract Module	13
Extract - Proc Module	15
Extract - ROI Module	18
Extract - Track Module	18
moseq2_extract.helpers package	20
Helpers - Data Module	20
Helpers - Extract Module	22
Helpers - Wrappers Module	23
moseq2_extract.io package	24
IO - Image Module	24
IO - Video Module	24
Index	27
Index	29
Python Module Index	35

Welcome to moseq2-extract's documentation!

moseq2_extract package

CLI Module

moseq2-extract

```
moseq2-extract [OPTIONS] COMMAND [ARGS]...
```

aggregate-results

Copies all extracted results (h5, yaml, avi) files from all extracted sessions to a new directory,

moseq2-extract aggregate-results [OPTIONS]

Options

-i, --input-dir <input_dir>
Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-extract/docs]

-f, --format <format>

New file name formats from resepective metadata [default: {start_time}_{session_name}_{subject_name}]

-o, --output-dir <output_dir>

Location for storing all results together [default: /Users/aymanzeine/Desktop/moseq/moseq2-extract/docs/aggregate_results/]

convert-raw-to-avi

Converts/Compresses a raw depth file into an avi file (with depth values) that is 8x smaller.

moseq2-extract convert-raw-to-avi [OPTIONS] INPUT_FILE

Options

-t, --threads <threads>

Number of threads for encoding [default: 3]

--delete

Delete raw file if encoding is sucessful [default: False]

--fps <fps>

Video FPS [default: 30]

-b, --chunk-size <chunk_size>

Chunk size [default: 3000]

-o, --output-file <output_file>

Path to output file

-v, --verbose <verbose>

Verbosity level out batch encoding. [0-1] [default: 0]

Arguments

INPUT_FILE

Required argument

copy-slice

Copies a segment of an input depth recording into a new video file.

moseq2-extract copy-slice [OPTIONS] INPUT_FILE

Options

-t, --threads <threads>

Number of threads for encoding [default: 3]

--delete

Delete raw file if encoding is sucessful [default: False]

--fps <fps>

Video FPS [default: 30]

-b, --chunk-size <chunk_size>

Chunk size [default: 3000]

-o, --output-file <output_file>

Path to output file

-c, --copy-slice <copy_slice>

Slice to copy [default: 0, 1000]

Arguments

INPUT_FILE

Required argument

download-flip-file

Downloads Flip-correction model that helps with orienting the mouse during extraction.

moseq2-extract download-flip-file [OPTIONS] [CONFIG_FILE]

Options

--output-dir <output_dir>

Temp storage [default: /Users/aymanzeine/moseq2]

Arguments

CONFIG_FILE

Optional argument

extract

Processes raw input depth recordings to output a cropped and oriented video of the mouse and saves the output+metadata to h5 files in the given output directory.

moseq2-extract extract [OPTIONS] INPUT_FILE

Options

- --config-file <config_file>
- --use-plane-bground

Use a plane fit for the background. Useful for mice that don't move much [default: False]

--output-dir <output_dir>

Output directory to save the results h5 file

--dilate-iterations <dilate_iterations>

Number of dilation iterations to increase bucket floor size. [default: 1]

- --bg-sort-roi-by-position-max-rois <bg_sort_roi_by_position_max_rois>
 Max original ROIs to sort by position [default: 2]
- --bg-sort-roi-by-position <bg_sort_roi_by_position>
 Sort ROIs by position [default: False]
- --bg-roi-fill-holes <bg_roi_fill_holes>
 Fill holes in ROI [default: True]
- --bg-roi-gradient-kernel <bg_roi_gradient_kernel>
 Kernel size for Sobel gradient filtering [default: 7]

```
--bg-roi-gradient-threshold <bg_roi_gradient_threshold>
  Gradient must be < this to include points [default: 3000]
--bg-roi-gradient-filter <bg_roi_gradient_filter>
  Exclude walls with gradient filtering [default: False]
--bg-roi-depth-range <bg_roi_depth_range>
  Range to search for floor of arena (in mm) [default: 650, 750]
--bg-roi-weights <bg_roi_weights>
  ROI feature weighting (area, extent, dist) [default: 1, 0.1, 1]
--bg-roi-index <bg_roi_index>
  Index of which background mask(s) to use [default: 0]
--bg-roi-shape <bg_roi_shape>
  Shape to use to dilate roi (ellipse or rect) [default: ellipse]
--bg-roi-dilate <bg_roi_dilate>
  Size of strel to dilate roi [default: 10, 10]
-c, --crop-size <crop_size>
 Width and height of cropped mouse image [default: 80, 80]
--min-height <min_height>
  Min mouse height from floor (mm) [default: 10]
--max-height <max_height>
  Max mouse height from floor (mm) [default: 100]
--detected-true-depth <detected_true_depth>
  Option to override automatic depth estimation during extraction. Either "auto" or a int value. [default: auto]
--fps <fps>
  Frame rate of camera [default: 30]
--flip-classifier <flip_classifier>
  Location of the flip classifier used to properly orient the mouse (.pkl file)
--flip-classifier-smoothing <flip_classifier_smoothing>
  Number of frames to smooth flip classifier probabilities [default: 51]
--use-tracking-model <use_tracking_model>
  Use an expectation-maximization style model to aid mouse tracking. Useful for data with cables [default: False]
--tracking-model-ll-threshold <tracking_model_ll_threshold>
  Threshold on log-likelihood for pixels to use for update during tracking [default: -100]
--tracking-model-mask-threshold <tracking_model_mask_threshold>
  Threshold on log-likelihood to include pixels for centroid and angle calculation [default: -16]
--tracking-model-ll-clip <tracking_model_ll_clip>
  Clip log-likelihoods below this value [default: -100]
--tracking-model-segment <tracking_model_segment>
  Segment likelihood mask from tracking model [default: True]
--tracking-model-init <tracking_model_init>
  Method for tracking model initialization [default: raw]
--cable-filter-iters <cable_filter_iters>
  Number of cable filter iterations [default: 0]
--cable-filter-shape <cable_filter_shape>
  Cable filter shape (rectangle or ellipse) [default: rectangle]
--cable-filter-size <cable_filter_size>
  Cable filter size (in pixels) [default: 5, 5]
--tail-filter-iters <tail filter iters>
  Number of tail filter iterations [default: 1]
--tail-filter-size <tail_filter_size>
```

```
Tail filter size [default: 9, 9]
--tail-filter-shape <tail_filter_shape>
 Tail filter shape [default: ellipse]
-s, --spatial-filter-size <spatial_filter_size>
  Space prefilter kernel (median filter, must be odd) [default: 3]
-t, --temporal-filter-size <temporal_filter_size>
  Time prefilter kernel (median filter, must be odd) [default: 0]
--chunk-size <chunk_size>
 Number of frames for each processing iteration [default: 1000]
--chunk-overlap <chunk_overlap>
  Frames overlapped in each chunk. Useful for cable tracking [default: 0]
--write-movie <write_movie>
  Write a results output movie including an extracted mouse [default: True]
--frame-dtype <frame_dtype>
  Data type for processed frames [default: uint8]
         Options: uint8|uint16
--centroid-hampel-span <centroid_hampel_span>
 Hampel filter span [default: 0]
--centroid-hampel-sig <centroid_hampel_sig>
  Hampel filter sig [default: 3]
--angle-hampel-span <angle_hampel_span>
  Angle filter span [default: 0]
--angle-hampel-sig <angle_hampel_sig>
  Angle filter sig [default: 3]
--model-smoothing-clips <model_smoothing_clips>
  Model smoothing clips [default: 0, 0]
--frame-trim <frame_trim>
  Frames to trim from beginning and end of data [default: 0, 0]
--compress <compress>
  Convert .dat to .avi after successful extraction [default: False]
--compress-chunk-size <compress_chunk_size>
  Chunk size for .avi compression [default: 3000]
--compress-threads <compress_threads>
 Number of threads for encoding [default: 3]
--verbose <verbose>
 Level of verbosity during extraction process. [0-2] [default: 0]
                                               Arguments
INPUT FILE
  Required argument
find-roi
Finds the ROI and background distance to subtract from frames when extracting.
```

```
moseq2-extract find-roi [OPTIONS] INPUT_FILE
```

Options

```
--config-file <config_file>
```

--use-plane-bground

Use a plane fit for the background. Useful for mice that don't move much [default: False]

- --output-dir <output_dir>
 Output directory to save the results h5 file
 --dilate-iterations <dilate_iterations>
 Number of dilation iterations to increase bucket floor size. [default: 1]
- --bg-sort-roi-by-position-max-rois <bg_sort_roi_by_position_max_rois>
 Max original ROIs to sort by position [default: 2]
- --bg-sort-roi-by-position <bg_sort_roi_by_position>
 Sort ROIs by position [default: False]
- --bg-roi-fill-holes <bg_roi_fill_holes>
 Fill holes in ROI [default: True]
- --bg-roi-gradient-kernel <bg_roi_gradient_kernel>
 Kernel size for Sobel gradient filtering [default: 7]
- --bg-roi-gradient-threshold <bg_roi_gradient_threshold>
 Gradient must be < this to include points [default: 3000]</pre>
- --bg-roi-gradient-filter <bg_roi_gradient_filter>
 Exclude walls with gradient filtering [default: False]
- --bg-roi-depth-range <bg_roi_depth_range>
 Range to search for floor of arena (in mm) [default: 650, 750]
- --bg-roi-weights <bg_roi_weights>
 ROI feature weighting (area, extent, dist) [default: 1, 0.1, 1]
- --bg-roi-index <bg_roi_index>
 Index of which background mask(s) to use [default: 0]
- --bg-roi-shape <bg_roi_shape>
 Shape to use to dilate roi (ellipse or rect) [default: ellipse]
- --bg-roi-dilate <bg_roi_dilate>
 Size of strel to dilate roi [default: 10, 10]

Arguments

INPUT_FILE

Required argument

generate-config

Generates a configuration file that holds editable options for extraction parameters.

moseq2-extract generate-config [OPTIONS]

Options

-o, --output-file <output_file>
[default: config.yaml]

generate-index

Generates an index YAML file containing all extracted session metadata.

moseq2-extract generate-index [OPTIONS]

Options

- -i, --input-dir <input_dir>
 - Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-extract/docs]
- -p, --pca-file <pca_file>
 - Path to PCA results [default: /Users/aymanzeine/Desktop/moseq/moseq2-extract/docs/ pca/pca scores.h5]
- -o, --output-file <output_file>

Location for storing index [default: /Users/aymanzeine/Desktop/moseq/moseq2-extract/docs/moseq2-index.yaml]

- -f, --filter <filter>
 Regex filter for metadata
- -a, --all-uuids <all_uuids>
 Use all uuids [default: False]
- --subpath <subpath>

Path substring to regulate paths included in an index file. [default: /proc/]

version

Print version number

moseq2-extract version [OPTIONS]

GUI Module

moseq2_extract.gui.aggregate_extract_results_command (input_dir, format, output_dir, subpath='/')
Finds all extracted h5, yaml and avi files and copies them all to a new directory relabeled with their respective session names. Also generates the index file.

Parameters:

- input_dir (str) (path to base directory to recursively search for h5s)
- format (str) (filename format for info to include in filenames)
- output_dir (str) (path to directory to save all aggregated results)

Returns: indexpath (str)

Return type: path to newly generated index file.

moseq2_extract.gui.check_progress (base_dir, progress_filepath)

Checks whether progress file exists and prompts user input on whether to overwrite, load old, or generate a new one

Parameters:

- base_dir (str) (path to directory to create/find progress file)
- progress_filepath (str) (path to progress filename)

Returns:

Return type: All restored variables or None.

moseq2_extract.gui.download_flip_command (output_dir, config_file=", selection=1)

Downloads flip classifier and saves its path in the inputted config file

Parameters:

- output_dir (str) (path to output directory to save flip classifier)
- config_file (str) (path to config file)
- selection (int) (index of which flip file to download (default is Adult male C57 classifer))

Returns:

Return type: None

moseq2_extract.gui.extract_command (input_file, output_dir, config_file, num_frames=None, skip=False)
Command to extract a full depth file

Parameters:

- input_file (str) (path to depthfile)
- output_dir (str) (path to output directory)
- config_file (str) (path to config file)
- num_frames (int) (number of frames to extract. All if None.)
- skip (bool) (skip already extracted file.)

Returns:

Return type: None

moseq2_extract.gui.extract_found_sessions (input_dir, config_file, ext, extract_all=True,
skip_extracted=False)

Searches for all depth files within input_directory with selected extension

Parameters:

- input_dir (str) (path to directory containing all session folders)
- config_file (str) (path to config file)
- ext (str) (file extension to search for)
- extract_all (bool) (if True, auto searches for all sessions, else, prompts user to select sessions individually.)
- skip extracted (bool) (indicates whether to skip already extracted session.)

Returns:

Return type: None

moseq2_extract.gui.find_roi_command (input_dir, config_file, exts=['dat', 'mkv', 'avi'], select_session=False,
default_session=0)

Computes ROI files given depth file. Will list out all available sessions to process and prompts user to input a corresponding session index to process.

Parameters:

- input_dir (str) (path to directory containing depth file)
- config_file (str) (path to config file)
- exts (list) (list of supported extensions)
- select_session (bool) (list all found sessions and allow user to select specific session to analyze via user-prompt)
- default_session (int) (index of the default session to find ROI for)

Returns: images (list of 2d arrays) (list of 2d array images to graph in Notebook.) filenames (list) (list of paths to respective image paths)

moseq2_extract.gui.generate_config_command (output_file)

Generates configuration file to use throughout pipeline.

Parameters: output_file (str) (path to saved config file.)

Returns: (str)

Return type: status message.

moseq2_extract.gui.generate_index_command (input_dir, pca_file, output_file, filter, all_uuids, subpath='/')
Generates Index File based on aggregated sessions

Parameters:

- input_dir (str) (path to aggregated_results/ dir)
- pca_file (str) (path to pca file)
- output_file (str) (index file name)
- filter (list) (keys to filter through)
- all uuids (list) (all extracted session uuids)

Returns: output_file (str)
Return type: path to index file.

moseq2_extract.gui.get_found_sessions (data_dir=", exts=['dat', 'mkv', 'avi'])
Find all depth recording sessions (with given extensions) to work on given base directory.

Parameters:

- data_dir (str) (path to directory containing all session folders)
- exts (list) (list of depth file extensions to search for)

Returns: data_dir (str) (path to base_dir to save in progress file) found_sessions (int) (number of found sessions with given extensions)

moseq2_extract.gui.restore_progress_vars (progress_file)

Restore all saved progress variables to Jupyter Notebook.

Parameters: progress_file (str) (path to progress file)

Returns:

Return type: All progress file variables

moseq2_extract.gui.sample_extract_command (input_dir, config_file, nframes, select_session=False,
default_session=0, exts=['dat', 'mkv', 'avi'])

Test extract command to extract a subset of the video.

Parameters:

- input_dir (str) (path to directory containing depth file to extract)
- config_file (str) (path to config file)
- nframes (int) (number of frames to extract)
- select_session (bool) (list all found sessions and allow user to select specific session to analyze via user-prompt)
- default_session (int) (index of the default session to find ROI for)
- exts (list) (list of supported depth file extensions.)

Returns: output_dir (str)

Return type: path to directory containing sample extraction results.

moseq2_extract.gui.update_progress (progress_file, varK, varV)
Updates progress file with new notebook variable

Parameters:

- progress_file (str) (path to progress file)
- varK (str) (key in progress file to update)
- varV (str) (updated value to write)

Returns:

Return type: None

moseq2_extract.gui.view_extraction (extractions, default=0)

Prompts user to select which extracted video(s) to preview.

Parameters:

- extractions (list) (list of paths to all extracted avi videos.)
- default (int) (index of the default extraction to display)

Returns: extractions (list)

Return type: list of selected extractions.

General Utilities Module

moseq2_extract.util.build_path (keys: dict, format_string: str, snake_case=True) → str

Produce a new file name using keys collected from extraction h5 files. The format string must be using python's formatting specification, i.e. '{subject_name}_{session_name}'.

Parameters:

- keys (dict) (dictionary specifying which keys used to produce the new file name)
- format_string (str) (the string to reformat using the keys dictionary)
- snake_case (bool) (whether to save the files with snake case)

Returns: out (str)

Return type: a newly formatted filename useable with any operating system

moseq2_extract.util.camel_to_snake(s)

Converts CamelCase to snake_case

Parameters: s (str) (CamelCase string to convert to snake_case.)

Returns: (str)

Return type: string in snake case

moseq2_extract.util.clean_dict (dct)

Standardizes types of dict value.

Parameters: dct (dict) (dict object with mixed type value objects.)

Returns: dct (dict)

Return type: dict object with list value objects.

moseq2_extract.util.clean_file_str (file_str: str, replace_with: str = '-') \rightarrow str Removes invalid characters for a file name from a string.

Parameters:

• file_str (str) (filename substring to replace)

• replace_with (str) (value to replace str with)

Returns: out (str)

Return type: cleaned file string

moseq2_extract.util.click_param_annot (click_cmd)

Given a click.Command instance, return a dict that maps option names to help strings. Currently skips click.Arguments, as they do not have help strings.

Parameters: click_cmd (click.Command) (command to introspect)

Returns: annotations (dict)

Return type: click.Option.human_readable_name as keys; click.Option.help as values

moseq2_extract.util.command_with_config (config_file_param_name)

moseq2_extract.util.convert_pxs_to_mm (coords, resolution=512, 424, field_of_view=70.6, 60, true_depth=673.1)

Converts x, y coordinates in pixel space to mm.

http://stackoverflow.com/questions/17832238/kinect-intrinsic-parameters-from-field-of-view/18199938#18199938 http://www.imaginativeuniversal.com/blog/post/2014/03/05/quick-reference-kinect-1-vs-kinect-2.aspx http://smeenk.com/kinect-field-of-view-comparison/

Parameters:

- coords (list) (list of x,y pixel coordinates)
- resolution (tuple) (image dimensions)
- field_of_view (tuple) (width and height scaling params)
- true depth (float) (detected true depth)

Returns: new_coords (list)

Return type: x,y coordinates in mm

moseq2_extract.util.convert_raw_to_avi_function (input_file, chunk_size=2000, fps=30, delete=False, threads=3)

Converts depth file to avi file.

Parameters:

- input file (str) (path to depth file)
- chunk_size (int) (size of chunks to process at a time)
- fps (int) (frames per second)
- delete (bool) (whether to delete original depth file)
- threads (int) (number of threads to write video.)

Returns:

Return type: None

moseq2_extract.util.dict_to_h5 (h5, dic, root='/', annotations=None)

Save an dict to an h5 file, mounting at root. Keys are mapped to group names recursively.

- h5 (h5py.File instance) (h5py.file object to operate on)
- dic (dict) (dictionary of data to write)
- root (string) (group on which to add additional groups and datasets)
- annotations (dict) (annotation data to add to corresponding h5 datasets. Should contain same keys as dic.)

Returns:

Return type: None

moseq2_extract.util.escape_path (path)

Given current path, will return a path to return to original base directory. (Used in recursive h5 search, etc.)

Parameters: path (str) (path to current working dir)

Returns: path (str)

Return type: path to original base_dir

moseq2_extract.util.gen_batch_sequence (nframes, chunk_size, overlap, offset=0) Generates batches used to chunk videos prior to extraction.

Parameters:

- nframes (int) (total number of frames)
- chunk_size (int) (desired chunk size)
- overlap (int) (number of overlapping frames)
- offset (int) (frame offset)

Returns:

Return type: Yields list of batches

moseq2_extract.util.get_bucket_center (img, true_depth, threshold=650)

https://stackoverflow.com/questions/19768508/python-opencv-finding-circle-sun-coordinates-of-center-the-circle-from-pictu Finds Centroid coordinates of circular bucket.

Parameters:

- img (2d np.ndaarray) (original background image.)
- true_depth (float) (distance value from camera to bucket floor (automatically pre-computed))
- threshold (float) (distance values to accept region into detected circle. (used to reduce fall noise interference))

Returns: cX (int) (x-coordinate of circle centroid) cY (int) (y-coordinate of circle centroid)

moseq2_extract.util.graduate_dilated_wall_area (bground_im, config_data, strel_dilate, true_depth,
output_dir)

Creates a gradient to represent the dilated (now visible) bucket wall regions. Only is used if background is dilated to capture larger rodents in convex shaped buckets (_/). This is done to handle noise attributed by bucket walls being slanted, and thus being picked up as large noise depth values. Moreover, to appropriately subtract the background from input images during extraction without obscuring the rodent, or including unwanted wall regions. :Parameters: * bground_im (2d np.ndarray) (the bucket floor image computed as the median distance throughout the session.)

- config_data (dict) (dictionary containing helper user configuration parameters.)
- strel_dilate (cv2.structuringElement) (dilation structuring element used to dilate background image.)
- true depth (float) (median distance computed throughout recording.)
- output_dir (str) (path to save newly computed background to use.)

Returns: bground im (2d np.ndarray)

Return type: the new background image with a gradient around the floor from high to low depth values.

 $moseq2_extract.util.h5_to_dict (h5file, path) \rightarrow dict$ Loads h5 contents to dictionary object.

• h5file (str or h5py.File) (file path to the given h5 file or the h5 file handle)

• path (str) (path to the base dataset within the h5 file)

Returns: out (dict)

Return type: a dict with h5 file contents with the same path structure

moseq2_extract.util.load_metadata (metadata_file)

Loads metadata.

Parameters: metadata_file (str) (path to metadata file)

Returns:

Return type: metadata (dict)

moseq2_extract.util.load_textdata (data_file, dtype=<class 'numpy.float32'>)
 Loads timestamp from txt/csv file

Parameters:

• data_file (str) (path to timestamp file)

dtype (dtype) (data type of timestamps)

Returns: data (np.ndarray) (timestamp data) timestamps (np.array) (time stamp keynames.)

moseq2_extract.util.load_timestamps (timestamp_file, col=0)

Read timestamps from space delimited text file.

Parameters:

• timestamp_file (str) (path to timestamp file)

• col (int) (column in ts file read.)

Returns: ts (list)

Return type: list of timestamps

moseq2_extract.util.make_gradient (width, height, h, k, a, b, theta=0)

https://stackoverflow.com/questions/49829783/draw-a-gradual-change-ellipse-in-skimage/49848093#49848093 Creates gradient around bucket floor representing slanted wall values. This is done by drawing an "ellipse" of equal x,y radii, resulting in a circle with weighted depth values from highest to lowest surrounding the circumference of the circle

Parameters:

• width (int) (bounding box width)

· height (int) bounding box height

• h (int) (centroid x coordinate)

• k (int) (centroid y coordinate)

• a (int) (x-radius of drawn ellipse)

• b (int) (y-radius of drawn ellipse)

• theta (float) (degree to rotate ellipse in radians. (has no effect if drawing a circle))

Returns:

(2d np.ndarray) (numpy array with weighted values from 0.08 -> 0.8 representing the proportion of values) to create a gradient from. 0.8 being the proportioned values closest to the circle wall.

moseg2 extract.util.mouse threshold filter (h5file, thresh=0)

Filters frames in h5 files by threshold value

Parameters:

• h5file (str) (path to h5 file)

• thresh (int) (threshold at which to apply filter)

Returns: (3d-np boolean array)

Return type: array of regions to include after threshold filter.

moseq2_extract.util.read_yaml (yaml_file)

Reads yaml file into dict object

Parameters: yaml_file (str) (path to yaml file)

Returns: return_dict (dict)
Return type: dict of yaml contents

moseq2_extract.util.recursive_find_h5s

(root_dir='/Users/aymanzeine/Desktop/moseq/moseq2-extract/docs', ext='.h5', yaml_string='{}.yaml')

Recursively find h5 files, along with yaml files with the same basename

Parameters:

- root_dir (str) (path to base directory to begin recursive search in.)
- ext (str) (extension to search for)

• yaml_string (str) (string for filename formatting when saving data)

Returns: h5s (list) (list of found h5 files) dicts (list) (list of found metadata files) yamls (list) (list of

found yaml files)

moseq2_extract.util.recursive_find_unextracted_dirs

(root_dir='/Users/aymanzeine/Desktop/moseq/moseq2-extract/docs', session_pattern='session_\\d+\\.(?:tgz|tar\\.gz)', filename='.dat', yaml_path='proc/results_00.yaml', metadata_path='metadata.json', skip_checks=True)

Recursively find unextracted (or incompletely extracted) directories

Parameters:

- root_dir (os Path-like) (path to base directory to start recursive search from.)
- session_pattern (str) (folder name pattern to search for)
- filename (str) (file extension to search for)
- yaml_path (str) (path to respective extracted metadata)
- metadata_path (str) (path to relative metadata.json files)
- **skip_checks (bool)** (indicates whether to check if the files exist at the given relative paths)

Returns: proc_dirs (1d-list)

Return type: list of paths to each unextracted session's proc/ directory

moseq2_extract.util.scalar_attributes ()

Gets scalar attributes

Returns: attributes (dict)

Return type: collection of metadata keys and descriptions.

moseq2_extract.util.select_strel (string='e', size=10, 10)

Returns structuring element of specified shape.

Parameters:

- string (str) (indicates whether to use ellipse or rectangle)
- size (tuple) (size of structuring element)

Returns:

Return type: strel (cv2.StructuringElement)

```
moseq2_extract.util.strided_app (a, L, S)
```

from https://stackoverflow.com/questions/40084931/taking-subarrays-from-numpy-array-with-given-stride-stepsize /40085052#40085052

Parameters:

- a (np.ndarray) array to get subarrarys from.
- L (int) Window Length
- S (int) Stride size

Returns:

Return type: (np.ndarray) - array of subarrays at stride S.

moseq2 extract.util.time str for filename (time str: str) → str

Process the time string supplied by moseq to be used in a filename. This removes colons, milliseconds, and timezones.

Parameters: time_str (str) (time str to format)

Returns: out (str)

Return type: formatted timestamp str

Subpackages

moseq2_extract.extract package

Extract - Extract Module

This function looks for a mouse in background-subtracted frames from a chunk of depth video. It is called from the moseq2_extract.helpers.extract module.

- chunk (3d np.ndarray) (chunk to extract (chunksize, height, width))
- use_em_tracker (bool) (The EM tracker uses expectation-maximization to fit a 3D gaussian on a frame-by-frame) basis to the mouse's body and determine if pixels are mouse vs cable.
- prefilter_space (tuple) (spatial kernel size)
- prefilter_time (tuple) (temporal kernel size)
- iters_tail (int) (number of filtering iterations on mouse tail)
- iters_min (int) (minimum tail filtering filter kernel size)
- strel_tail (cv2::StructuringElement Ellipse) (filtering kernel size to filter out mouse tail.)
- strel_min (cv2::StructuringElement Rectangle) (filtering kernel size to filter mouse body in cable recording cases.)
- min_height (int) (minimum (mm) distance of mouse to floor.)
- max height (int) (maximum (mm) distance of mouse to floor.)
- mask_threshold (int) (Threshold on log-likelihood to include pixels for centroid and angle calculation)
- use_cc (bool) (boolean to use connected components in cv2 structuring elements)
- bground (np.ndarray) (numpy array represented previously computed background)
- roi (np.ndarray) (numpy array represented previously computed roi)
- rho_mean (int) (smoothing parameter for the mean)
- rho_cov (int) (smoothing parameter for the covariance)
- tracking_II_threshold (float) (threshold for calling pixels a cable vs a mouse (usually between -16 to -12).) – If the log-likelihood falls below this value, pixels are considered cable.
- tracking_segment (bool) (boolean for whether to use only the largest blob for EM updates.)
- tracking_init_mean (float) (Initialized mean value for EM Tracking)
- tracking_init_cov (float) (Initialized covariance value for EM Tracking)
- tracking_init_strel (cv2::StructuringElement Ellipse)
- flip_classifier (str) (path to pre-selected flip classifier.)
- flip_smoothing (int) (amount of smoothing to use for flip classifier.)
- frame dtype (str) (Data type for processed frames)
- save path ((str): Path to save extracted results)
- progress_bar (bool) (Display progress bar)
- crop_size (tuple) (size of the cropped mouse image.)
- true_depth (float) (previously computed detected true depth value.)
- centroid_hampel_span (int) (Hampel filter span kernel size)
- centroid_hampel_sig (int) (Hampel filter standard deviation)
- angle_hampel_span (int) (Angle filter span kernel size)
- angle_hampel_sig (int) (Angle filter standard deviation)
- model_smoothing_clips (tuple) (Model smoothing clips)
- tracking_model_init (str) (Method for tracking model initialization)
- verbose (int) (Level of verbosity during extraction process. [0-2])

Returns: results ((3d np.ndarray) - (nframes, crop height, crop width)) extracted cropped, oriented

and centered RGB video chunk to be written to file.

Extract - Proc Module

moseq2_extract.extract.proc.apply_roi (frames, roi)
Apply ROI to data, consider adding constraints (e.g. mod32==0).

Parameters:

• frames (3d np.ndarray) (input frames to apply ROI.)

• roi (2d np.ndarray) (selected ROI to extract from input images.)

Returns: cropped_frames (3d np.ndarray)

Return type: Frames cropped around ROI Bounding Box.

moseq2_extract.extract.proc.clean_frames (frames, prefilter_space=3, prefilter_time=None, strel_tail=array([[0, 0, 0, 1, 0, 0, 0], [0, 1, 1, 1, 1, 0], [1, 1, 1, 1, 1, 1], [1, 1, 1, 1, 1, 1], [1, 1, 1, 1, 1], [1, 1, 1, 1, 1], [0, 1, 1, 1, 1, 0], [0, 0, 0, 1, 0, 0, 0]], dtype=uint8), iters_tail=None, frame_dtype='uint8', strel_min=array([[1, 1, 1, 1, 1], [1, 1, 1, 1, 1], [1, 1, 1, 1, 1], [1, 1, 1, 1, 1]], dtype=uint8), iters_min=None, progress_bar=True, qui=False, verbose=0)

Simple filtering, median filter and morphological opening.

Parameters:

- frames (3d np.ndarray) (Frames (nframes x r x c) to filter.)
- prefilter_space (tuple) (kernel size for spatial filtering)
- prefilter_time (tuple) (kernel size for temporal filtering)
- strel_tail (cv2.StructuringElement) (Element for tail filtering.)
- iters_tail (int) (number of iterations to run opening)
- frame_dtype (str) (frame encodings)
- strel_min (int) (minimum kernel size)
- iters_min (int) (minimum number of filtering iterations)
- progress_bar (bool) (display progress bar)
- gui (bool) (indicate GUI is executing function)
- verbose (bool) (display progress)

Returns: filtered_frames (3d np array)

Return type: frame x r x c

moseq2_extract.extract.proc.compute_scalars (frames, track_features, min_height=10, max_height=100, true_depth=673.1)

Computes scalars.

Parameters:

- frames (3d np.ndarray) (frames x r x c, uncropped mouse)
- track features (dict) (dictionary with tracking variables (centroid and orientation))
- min_height (float) (minimum height of the mouse)
- max_height (float) (maximum height of the mouse)
- true_depth (float) (detected true depth)

Returns: features (dict)

Return type: dictionary of scalars

moseq2_extract.extract.proc.crop_and_rotate_frames (frames, features, crop_size=80, 80, progress_bar=True, gui=False, verbose=0)

Crops mouse from image and orients it s.t it is always facing east.

- frames (3d np.ndarray) (frames to crop and rotate)
- features (dict) (dict of extracted features, found in result_00.h5 files.)
- crop_size (tuple) (size of cropped image.)
- progress_bar (bool) (Display progress bar.)
- gui (bool) (indicate GUI is executing function)
- verbose (bool) (display progress)

Returns: cropped_frames (3d np.ndarray)

Return type: Crop and rotated frames.

moseq2_extract.extract.proc.feature_hampel_filter (features, centroid_hampel_span=None, centroid_hampel_sig=3, angle_hampel_span=None, angle_hampel_sig=3)

Filters computed extraction features using Hampel Filtering.

Parameters:

- features (dict) (dictionary of video features)
- centroid_hampel_span (int) (Centroid Hampel Span Filtering Kernel Size)
- centroid_hampel_sig (int) (Centroid Hampel Signal Filtering Kernel Size)
- angle_hampel_span (int) (Angle Hampel Span Filtering Kernel Size)
- angle_hampel_sig (int) (Angle Hampel Span Filtering Kernel Size)

Returns: features (dict)

Return type: filtered version of input dict.

moseq2_extract.extract.proc.get_bbox (roi)

Given a binary mask, return an array with the x and y boundaries

Parameters: roi (2d np.ndarray) (ROI boolean mask to calculate bounding box.)

Returns: bbox (2d np.ndarray)
Return type: Bounding Box around ROI

moseq2_extract.extract.proc.get_bground_im (frames)

Returns background

Parameters: frames (3d numpy array) (frames x r x c, uncropped mouse)

Returns: bground (2d numpy array)

Return type: r x c, background image

moseq2_extract.extract.proc.get_bground_im_file (frames_file, frame_stride=500, med_scale=5,
**kwargs)

Returns background from file

Parameters:

- frames_file (str) (path to data with frames)
- frame_stride (int) (stride size between frames for median bground calculation)
- med_scale (int) (kernel size for median blur for background images.)
- kwargs

Returns: bground (2d numpy array)
Return type: r x c, background image

logged? extract extract proc get fling (frames flin file-None

moseq2_extract.extract.proc.get_flips (frames, flip_file=None, smoothing=None)
Predicts frames where mouse orientation is flipped to later correct.

Parameters:

- frames (3d numpy array) (frames x r x c, cropped mouse)
- flip_file (str) (path to joblib dump of scipy random forest classifier)
- smoothing (int) (kernel size for median filter smoothing of random forest probabilities)

Returns: flips (bool array)

Return type: true for flips

moseq2_extract.proc.get_frame_features (frames, frame_threshold=10, mask=array([], dtype=float64), mask_threshold=- 30, use_cc=False, progress_bar=True, gui=False, verbose=0)

Use image moments to compute features of the largest object in the frame

Parameters:

- frames (3d np.ndarray) (input frames)
- frame_threshold (int) (threshold in mm separating floor from mouse)
- mask (3d np.ndarray) (input frame mask for parts not to filter.)
- mask_threshold (int) (threshold to include regions into mask.)
- use_cc (bool) (Use connected components.)
- progress_bar (bool) (Display progress bar.)
- gui (bool) (indicate GUI is executing function)
- verbose (bool) (display progress)

Returns: features (dict of lists) (dictionary with simple image features) mask (3d np.ndarray) (input frame mask.)

moseq2_extract.extract.proc.get_largest_cc (frames, progress_bar=False)
Returns largest connected component blob in image

Parameters:

- frames (3d numpy array) (frames x r x c, uncropped mouse)
- progress_bar (bool) (display progress bar)

Returns: flips (3d bool array)

Return type: frames x r x c, true where blob was found

Parameters:

- depth_image (2d np.ndarray) (Singular depth image frame.)
- strel_dilate (cv2.StructuringElement Rectangle) (dilation shape to use.)
- dilate iters (int) (number of dilation iterations.)
- strel_erode (int) (image erosion kernel size.)
- noise_tolerance (int) (threshold to use for noise filtering.)
- weights (tuple) (weights describing threshold to accept ROI.)
- overlap roi (np.ndarray) (list of ROI boolean arrays to possibly combine.)
- gradient_filter (bool) (Boolean for whether to use a gradient filter.)
- gradient_kernel (tuple) (Kernel size of length 2, e.g. (1, 1.5))
- gradient_threshold (int) (Threshold for noise gradient filtering)
- fill_holes (bool) (Boolean to fill any missing regions within the ROI.)
- gui (bool) (Boolean for whether function is running on GUI.)
- verbose (bool) (Boolean for whether to display progress)
- kwargs

Returns: rois (list) (list of 2d roi images.) roi plane (2d np.ndarray) (computed ROI Plane using

RANSAC.) bboxes (list) (list of computed bounding boxes for each respective ROI.) label_im (list) (list of scikit-image image properties) ranks (list) (list of ROI ranks.)

shape_index (list) (list of rank means.)

moseq2_extract.extract.proc.im_moment_features (IM)

Use the method of moments and centralized moments to get image properties.

Parameters: IM (2d numpy array) (depth image)

Returns: features (dict) – centroid, and ellipse axis length

Return type: returns a dictionary with orientation,

moseq2_extract.extract.proc.model_smoother (features, II=None, clips=- 300, - 125) Spatial feature filtering.

Parameters:

• features (dict) (dictionary of extraction scalar features)

• II (np.array) (list of loglikelihoods of pixels in frame)

• clips (tuple) (tuple to ensure video is indexed properly)

Returns:

Return type: features (dict) - smoothed version of input features

Extract - ROI Module

moseq2_extract.extract.roi.plane_fit3 (points)

Fit a plane to 3 points (min number of points for fitting a plane)

Parameters: points (2d numpy array) (each row is a group of points, columns correspond to x,y,z.)

Returns: plane (1d numpy array)

Return type: linear plane fit—>a*x+b*y+c*z+d

moseq2_extract.extract.roi.plane_ransac (depth_image, depth_range=650, 750, iters=1000, noise_tolerance=30, in_ratio=0.1, progress_bar=True, mask=None, gui=False, verbose=0)

Naive RANSAC implementation for plane fitting

Parameters:

- depth_image (2d numpy array) (hxw, background image to fit plane to)
- depth_range (tuple) (min/max depth (mm) to consider pixels for plane)
- iters (int) (number of RANSAC iterations)
- noise_tolerance (float) (dist. from plane to consider a point an inlier)
- in_ratio (float) (frac. of points required to consider a plane fit good)
- progress_bar (bool) (display progress bar)
- mask (bool 2d np.array) (boolean mask to find region to use)
- gui (bool) (whether GUI is used.)
- verbose (int) (0 or 1; 1 to print all information.)

Returns: best_plane (1d numpy array) (plane fit to data) dist (1d numpy array) (distance of the

calculated coordinates and "best plane")

Extract - Track Module

moseq2_extract.extract.track.em_get_11 (frames, mean, cov, progress_bar=True)
Returns likelihoods for each frame given tracker parameters

- frames (3d numpy array) (depth frames)
- mean (2d numpy array) (frames x d, mean estimates)
- cov (3d numpy array) (frames x d x d, covariance estimates)
- progress_bar (bool) (use a progress bar)

Returns: II (3d numpy array)

Return type: frames x rows x columns, log likelihood of all pixels in each frame

Parameters:

- depth_frame (2d numpy array) (depth frame to initialize mask with.)
- depth_floor (float) (distance from camera to bucket floor.)
- depth_ceiling (float) (max depth value.)
- init_strel (cv2.structuringElement) (structuring Element to compute mask.)
- strel_iters (int) (number of morphological iterations.)

Returns: mouse_mask (2d numpy array)

Return type: mask of depth frame.

moseq2_extract.extract.track.em_iter (data, mean, cov, lamd=0.1, epsilon=0.1, max_iter=25)
Single iteration of EM tracker

Parameters:

- data (3d numpy array) (nx3, x, y, z coordinates to use)
- mean (1d numpy array) (dx1, current mean estimate)
- cov (2d numpy array) (dxd, current covariance estimate)
- lambd (float) (constant to add to diagonal of covariance matrix)
- epsilon (float) (tolerance on change in likelihood to terminate iteration)
- max_iter (int) (maximum number of EM iterations)

Returns: mean (1d numpy array) (updated mean) cov (2d numpy array) (updated covariance)

Naive tracker, use EM update rules to follow a 3D Gaussian

around the room.

- frames (3d numpy array) (filtered frames nframes x r x c.)
- raw_frames (3d numpy array) (chunk to track mouse in.)
- segment (bool) (use only the largest blob for em updates)
- II_threshold (float) (threshold on log likelihood for segmentation)
- rho_mean (float) (smoothing parameter for the mean)
- rho_cov (float) (smoothing parameter for the covariance)
- depth_floor (float) (height in mm for separating mouse from floor)
- depth_ceiling (float) (max height in mm for mouse from floor.)
- progress_bar (bool) (display progress bar.)
- init_mean (np.ndarray) (array of inital frame pixel means.)
- init_cov (np.ndarray) (array of inital frame pixel covariances.)
- init_frames (int) (number of frames to include in the init calulation)
- init_method (str) (mode in which to process inputs)
- init_strel (cv2.structuringElement) (structuring Element to compute mask.)

Returns: model_parameters (dict)

Return type: mean and covariance estimates for each frame

moseq2_extract.helpers package

Helpers - Data Module

moseq2_extract.helpers.data.build_manifest (loaded, format, snake_case=True) aggregate_results() Helper Function. Builds a manifest file used to contain extraction result metadata from h5 and yaml files.

Parameters:

- loaded (list of dicts) (list of dicts containing loaded h5 data.)
- format (str) (filename format indicating the new name for the metadata files in the aggregate_results dir.)
- snake case (bool) (whether to save the files using snake case)

Returns: manifest (dict)

Return type: dictionary of extraction metadata.

moseq2_extract.helpers.data.copy_manifest_results (manifest, output_dir)
Copies all considated manifest results to their respective output files.

Parameters:

- manifest (dict) (manifest dictionary containing all extraction h5 metadata to save)
- output_dir (str) (path to directory where extraction results will be aggregated.)

Returns:

Return type: None

moseq2_extract.helpers.data.create_extract_h5 (f, acquisition_metadata, config_data, status_dict, scalars, scalars_attrs, nframes, true_depth, roi, bground_im, first_frame, timestamps, extract=None)

Creates h5 file that holds all extracted frames and other metadata (such as scalars).

- f (h5py.File object) (opened h5 file object to write to.)
- acquisition_metadata (dict) (Dictionary containing extracted session acquisition metadata.)
- config_data (dict) (dictionary object containing all required extraction parameters. (auto generated))
- status_dict (dict) (dictionary that helps indicate if the session has been extracted fully.)
- scalars (list) (list of computed scalar metadata.)
- scalars_attrs (dict) (dict of respective computed scalar attributes and descriptions to save.)
- nframes (int) (number of frames being recorded)
- true_depth (float) (computed detected true depth)
- roi (2d np.ndarray) (Computed 2D ROI Image.)
- bground_im (2d np.ndarray) (Computed 2D Background Image.)
- first_frame (2d np.ndarray) (Computed 2D First Frame Image.)
- timestamps (np.array) (Array of session timestamps.)
- extract (moseq2_extract.cli.extract function) (Used to preseve CLI state parameters in extraction h5.)

Returns:

Return type: None

moseq2_extract.helpers.data.get_selected_sessions (to_extract, extract_all)

Given user input, the function will return either selected sessions to extract, or all the sessions.

Parameters:

- to_extract (list) (list of paths to sessions to extract)
- extract_all (bool) (boolean to include all sessions and skip user-input prompt.)

Returns: to_extract (list)

Return type: new list of selected sessions to extract.

moseq2_extract.helpers.data.handle_extract_metadata (input_file, dirname, config_data, nframes) Extracts metadata from input depth files, either raw or compressed.

Parameters:

- input_file (str) (path to input file to extract)
- dirname (str) (path to directory where extraction files reside.)
- config_data (dict) (dictionary object containing all required extraction parameters. (auto generated))
- **nframes (int)** (number of frames to extract.)

Returns:

metadata_path (str) (path to respective metadata.json) timestamp_path (str) (path to respective depth_ts.txt or similar) alternate_correct (bool) (indicator for whether an alternate timestamp file was used) tar (bool) (indicator for whether the file is compressed.) nframes (int) (number of frames to extract) first_frame_idx (int) (index number of first frame in extraction.) last frame idx (int) (index number of last frame in extraction)

moseq2_extract.helpers.data.load_h5s (to_load, snake_case=True)
aggregate_results() Helper Function to load h5 files.

Parameters:

- to_load (list) (list of paths to h5 files.)
- snake case (bool) (whether to save the files using snake case)

Returns: loaded (list)

Return type: list of loaded h5 dicts.

Helpers - Extract Module

moseq2_extract.helpers.extract.process_extract_batches (f, input_file, config_data, bground_im, roi, scalars, frame_batches, first_frame_idx, true_depth, tar, strel_tail, strel_min, output_dir, output_filename)

Compute extracted frames and save them to h5 files and avi files.

Parameters:

- f (h5py.File) (opened h5 file to write extracted batches to)
- input_file (str) (path to depth file)
- config_data (dict) (dictionary containing extraction parameters (autogenerated))
- bground_im (2d numpy array) (r x c, background image)
- roi (2d numpy array) (r x c, roi image)
- scalars (list) (list of keys to scalar attribute values)
- frame_batches (list) (list of batches of frames to serially process.)
- first_frame_idx (int) (index of starting frame.)
- true_depth (float) (computed detected true depth.)
- tar (bool) (compressed file indicator.)
- strel_tail (cv2.StructuringElement) (Element for tail filtering.)
- strel_min (int) (minimum kernel size)
- output_dir (str) (path to output directory that contains the extracted data, e.g. (proc/).)
- output_filename (str) (name of h5 file containing extraction data, e.g. (results_00).)

Returns: video_pipe (bool)

Return type: boolean for whether function is done writing to video file.

moseq2_extract.helpers.extract.run_local_extract (to_extract, params, prefix, skip_extracted=False)
Runs the extract command on given list of sessions to extract on local platform. This function is meant for the GUI interface to utilize the moseq2-batch extract functionality.

Parameters:

- to_extract (list) (list of paths to files to extract)
- params (dict) (dictionary of ROI metadata from config file.)
- prefix (str) (prefix to CLI extraction command.)
- skip_extracted (bool) (Whether to skip already extracted session.)

Returns:

Return type: None

moseq2_extract.helpers.extract.run_slurm_extract (to_extract, params, partition, prefix, escape_path,
skip_extracted=False)

Runs the extract command on given list of sessions to extract on SLURM platform. This function is meant for the GUI interface to utilize the moseg2-batch extract functionality.

Parameters:

- to_extract (list) (list of paths to files to extract)
- params (dict) (dictionary of ROI metadata from config file.)
- partition (str) (name of slurm partition to use)
- prefix (str) (prefix to CLI extraction command.)
- escape_path (function) (gets path to return to original base directory)
- skip_extracted (bool) (Whether to skip already extracted session.)

Returns:

Return type: None

Helpers - Wrappers Module

moseq2_extract.helpers.wrappers.aggregate_extract_results_wrapper (input_dir, format,
output_dir)

Copies all the h5, yaml and avi files generated from all successful extractions to a new directory to hold all the necessary data to continue down the moseq pipeline.

Parameters:

- input dir (str) (path to base directory containing all session folders)
- format (str) (string format for metadata to use as the new aggregated filename)
- output dir (str) (name of the directory to create and store all results in)

Returns:

Return type: None

moseq2_extract.helpers.wrappers.copy_h5_metadata_to_yaml_wrapper (input_dir,
h5_metadata_path)

Copy's user specified metadata from h5path to a yaml file.

Parameters:

- input_dir (str) (path to directory containing h5 files)
- h5_metadata_path (str) (path within h5 to desired metadata to copy to yaml.)

Returns:

Return type: None

moseq2_extract.helpers.wrappers.extract_wrapper (input_file, output_dir, config_data, num_frames=None, skip=False, extract=None, gui=False)

Wrapper function to run extract function for both GUI and CLI.

Parameters:

- input_file (str) (path to depth file)
- output_dir (str) (path to directory to save results in.)
- config_data (dict) (dictionary containing extraction parameters.)
- num_frames (int) (number of frames to extract. All if None.)
- skip (bool) (indicates whether to skip file if already extracted)
- extract (function) (extraction function state (Only passed by CLI))
- gui (bool) (indicates if GUI is running.)

Returns: output_dir (str)

Return type: path to directory containing extraction (only if gui==True)

moseq2_extract.helpers.wrappers.flip_file_wrapper (config_file, output_dir, selected_flip=1,
qui=False)

Wrapper function to download and save flip classifiers. :Parameters: * config_file (str) (path to config file)

- output dir (str) (path to directory to save classifier in.)
- selected_flip (int) (index of desired flip classifier.)
- gui (bool) (indicates if the GUI is running.)

Returns:

Return type: None

moseq2_extract.helpers.wrappers.generate_index_wrapper (input_dir, pca_file, output_file, filter, all uuids, subpath='/proc/')

Generates index file containing a summary of all extracted sessions.

- input_dir (str) (directory to search for extracted sessions.)
- pca_file (str) (path to pca_scores file.)
- output_file (str) (preferred name of the index file.)
- filter (list) (list of metadata keys to conditionally filter.)
- all_uuids (list) (list of all session uuids.)

Returns: output_file (str)
Return type: path to index file.

moseq2_extract.helpers.wrappers.get_roi_wrapper (input_file, config_data, output_dir=None,
gui=False, extract_helper=False)

Wrapper function to compute ROI given depth file.

Parameters:

- input_file (str) (path to depth file.)
- config_data (dict) (dictionary of ROI extraction parameters.)
- output_dir (str) (path to desired directory to save results in.)
- gui (bool) (indicate whether GUI is running.)
- extract_helper (bool) (indicate whether this is being run independently or by extract function)

Returns:

if gui – output_dir (str): path to saved ROI results elif extract_helper – roi (2d array): ROI image to plot in GUI bground_im (2d array): Background image to plot in GUI first_frame (2d array): First frame image to plot in GUI

moseq2_extract.io package

10 - Image Module

moseq2_extract.io.image.read_image (filename, dtype='uint16', scale=True, scale_key='scale_factor')

Load image data, possibly with scale factor...

filename (str): path to file to write to.

image (2d numpy array): image to write scale (bool): indicates whether to scale image scale_key (str): indicates scale factor.

image (2d np array): loaded image

moseq2_extract.io.image.write_image (filename, image, scale=True, scale_factor=None, dtype='uint16',
metadata={}, compress=0)

Save image data, possibly with scale factor for easy display.

Parameters:

- filename (str) (path to file to write to.)
- image (2d numpy array) (the (unscaled) 2-D image to save)
- scale (bool) (flag to scale the image between the bounds of dtype)
- scale_factor (int) (factor by which to scale image)
- dtype (str) (array data type)
- metadata (dict) ([UNUSED] dictionary object that contains scaling info)
- compress (int) (image compression level)

Returns:

Return type: None

IO - Video Module

moseq2_extract.io.video.convert_mkv_to_avi (filename)

Converts Azure MKV video file format to AVI.

Parameters: filename (str) path to mkv file to convert

Returns: outpath (str)

Return type: path to converted AVI video file.

moseq2_extract.io.video.get_movie_info (filename, frame_dims=512, 424, bit_depth=16)
Returns dict of movie metadata.

Parameters:

• filename (str) (path to video file)

• frame_dims (tuple) (video dimensions)

bit_depth (int) (integer indicating data type encoding)

Returns: metadata (dict)

Return type: dictionary containing video file metadata

moseq2_extract.io.video.get_raw_info (filename, bit_depth=16, frame_dims=512, 424)
Gets info from a raw data file with specified frame dimensions and bit depth.

Parameters:

• filename (string) (name of raw data file)

• bit_depth (int) (bits per pixel (default: 16))

• frame_dims (tuple) (wxh or hxw of each frame)

Returns: file_info (dict)

Return type: dictionary containing depth file metadata

moseq2_extract.io.video.get_video_info (filename)

Get dimensions of data compressed using ffv1, along with duration via ffmpeg.

Parameters: filename (string) (name of file)

Returns: (dict)

Return type: dictionary containing video file metadata

moseq2_extract.io.video.load_movie_data (filename, frames=None, frame_dims=512, 424, bit_depth=16,
**kwargs)

Reads in frames

moseq2_extract.io.video.read_frames (filename, frames=range(0, 0), threads=6, fps=30, pixel_format='gray16le', frame_size=None, slices=24, slicecrc=1, get_cmd=False)

Reads in frames from the .nut/.avi file using a pipe from ffmpeg.

Parameters:

- filename (str) (filename to get frames from)
- frames (list or 1d numpy array) (list of frames to grab)
- threads (int) (number of threads to use for decode)
- fps (int) (frame rate of camera in Hz)
- pixel format (str) (ffmpeg pixel format of data)
- frame_size (str) (wxh frame size in pixels)
- slices (int) (number of slices to use for decode)
- slicecrc (int) (check integrity of slices)
- get_cmd (bool) (indicates whether function should return ffmpeg command (instead of executing).)

Returns: video (3d numpy array)

Return type: frames x h x w

moseq2_extract.io.video.read_frames_raw (filename, frames=None, frame_dims=512, 424, bit_depth=16, dtype='<i2', tar_object=None)</pre>

Reads in data from raw binary file.

- filename (string) (name of raw data file)
- frames (list or range) (frames to extract)
- frame_dims (tuple) (wxh of frames in pixels)
- bit_depth (int) (bits per pixel (default: 16))

• tar_object (tarfile.TarFile) (TarFile object, used for loading data directly from tgz)

Returns: chunk (numpy ndarray)

Return type: nframes x h x w

moseq2_extract.io.video.write_frames (filename, frames, threads=6, fps=30, pixel_format='gray16le', codec='ffv1', close_pipe=True, pipe=None, slices=24, slicecrc=1, frame_size=None, get_cmd=False, verbose=0) Write frames to avi file using the ffv1 lossless encoder

Parameters:

- filename (str) (path to file to write to.)
- frames (np.ndarray) (frames to write)
- threads (int) (number of threads to write video)
- fps (int) (frames per second)
- pixel_format (str) (format video color scheme)
- codec (str) (ffmpeg encoding-writer method to use)
- close_pipe (bool) (indicates to close the open pipe to video when done writing.)
- pipe (subProcess.Pipe) (pipe to currently open video file.)
- slices (int) (number of frame slices to write at a time.)
- slicecrc (int) (check integrity of slices)
- frame_size (tuple) (shape/dimensions of image.)
- get_cmd (bool) (indicates whether function should return ffmpeg command (instead of executing))
- verbose (bool) (output progress.)

Returns: pipe (subProcess.Pipe)

Return type: indicates whether video writing is complete.

moseq2_extract.io.video.write_frames_preview (filename, frames=array([], dtype=float64), threads=6, fps=30, pixel_format='rgb24', codec='h264', slices=24, slicecrc=1, frame_size=None, depth_min=0, depth_max=80, get_cmd=False, cmap='jet', pipe=None, close_pipe=True, frame_range=None)

Writes out a false-colored mp4 video.

- filename (str) (path to file to write to.)
- frames (np.ndarray) (frames to write)
- threads (int) (number of threads to write video)
- fps (int) (frames per second)
- pixel_format (str) (format video color scheme)
- codec (str) (ffmpeg encoding-writer method to use)
- slices (int) (number of frame slices to write at a time.)
- slicecrc (int) (check integrity of slices)
- frame_size (tuple) (shape/dimensions of image.)
- depth_min (int) (minimum mouse depth from floor in (mm))
- depth_max (int) (maximum mouse depth from floor in (mm))
- get_cmd (bool) (indicates whether function should return ffmpeg command (instead of executing))
- cmap (str) (color map to use.)
- pipe (subProcess.Pipe) (pipe to currently open video file.)
- close_pipe (bool) (indicates to close the open pipe to video when done writing.)
- frame_range (range()) (frame indices to write on video)

Returns: pipe (subProcess.Pipe)

Return type: indicates whether video writing is complete.

Index

• genindex

Index

Symbols

```
moseq2-extract-general hodexilter-iters <cable_filter_iters>
                                                                                                       moseq2-extract-ext
   --all-uuids <all_uuids>
                                                                                                       command line optic
                                  command line option
                                              moseq2catract-extract-shape <cable_filter_shape>
                                                                                                       moseq2-extract-ext
-angle-hampel-sig <angle_hampel_sig>
                                                                                                       command line optic
                                              command line option
                                               moseq2cable_filter_size <cable_filter_size>
                                                                                                      moseg2-extract-extra
angle-hampel-span <angle_hampel_span>
                                                                                                      command line option
                                               command line option
                                                moseq2texpact-expact-sig <centroid_hampel_sig>
                                                                                                        moseq2-extract-ex
g-roi-depth-range <bg_roi_depth_range>
                                                                                                        command line opt
                                                command line option
                                                --centroid-hampel-span <centroid_hampel_span>moseq2-extract-find-roi
                                                                                                         moseq2-extract-
                                                                                                         command line or
                                                command line option
                                          --chunk-overlap <chunk_overlap>
                                                                                                 moseq2-extract-extract
  --bg-roi-dilate <bg_roi_dilate>
                                                                                                 command line option
                                          command line option
                                                                                        moseq2-extract-convert-raw-to-avi
                                                          chunk-size <chunk size>
                                          moseq2-extract-find-roi
                                                                                        command line option
                                          command line option
                                                                                        moseq2-extract-copy-slice
                                               moseq2-extract-extract
bg-roi-fill-holes <bg_roi_fill_holes>
                                                                                        command line option
                                               command line option
                                                                                        moseq2-extract-extract command
                                               moseq2-extract-find-roi
                                                                                        line option
                                               command line option
                                                                                           moseq2-extract-extract
                                                                            <compress>
                                                     moseq2-extract-extract
i-gradient-filter <bg_roi_gradient_filter>
                                                                                           command line option
                                                     command line option
                                                                                                        moseq2-extract-ex
                                                    compress-chunk-size <compress_chunk_size>
                                                     moseq2-extract-find-roi
                                                                                                        command line opt
                                                     command line option
--compress-threads <compress_threads>
moseq2-extract-extract
                                                                                                    moseq2-extract-extract
i-gradient-kernel <bg_roi_gradient_kernel>
                                                                                                     command line option
                                                     command line option
    --config-file <config_file>
                                                                                              moseq2-extract-extract
                                                     moseq2-extract-find-roi
                                                                                              command line option
                                                     command line option
                                                                                              moseq2-extract-find-roi
                                                         moseq2-extract-extract
radient-threshold <bg_roi_gradient_threshold>
                                                                                              command line option
                                                         command line option
   --copy-slice <copy_slice>
                                                                                            moseq2-extract-copy-slice
                                                         moseq2-extract-find-roi
                                                                                            command line option
                                                         command line option <crop_size>
                                                                                            moseq2-extract-extract
                                                                                            command line option
                                         moseg2-extract-extract
   --bg-roi-index <bg_roi_index>
                                         command line option
                                                                                moseg2-extract-convert-raw-to-avi
                                                               --delete
                                                                                command line option
                                         moseq2-extract-find-roi
                                         command line option
                                                                                moseq2-extract-copy-slice
                                                                                command line option
                                         moseq2-extract-extract
   --bg-roi-shape <bg_roi_shape>
                                         command_line_optioned-true-depth <detected_true_depth>
                                                                                                        moseq2-extract-ex
                                                                                                        command line opt
                                         moseq2-extract-find-roi
                                                                                                     moseg2-extract-extra
                                         command line dpflarte-iterations <dilate_iterations>
                                                                                                      command line option
 --bg-roi-weights <bg_roi_weights>
                                           moseq2-extract-extract
                                           command line option
                                                                                                      moseq2-extract-find-
                                                                                                      command line option
                                           moseq2-extract-find-roi
                                                                                    moseq2-extract-generate-index
                                           command line option-filter <filter>
                                                                                    command line option
t-roi-by-position <bg_sort_roi_by_position>
                                                      moseq2-extract-extract
                                                      commignationer <flip_classifier>
                                                                                                   moseq2-extract-extract
                                                                                                   command line option
                                                      moseq2-extract-find-roi
                                                                                                               moseq2-e
                                          --flip-classifier-smoothing <flip_classifier_smoothing>
```

--format <format>

command

moseq2-extract-aggregate-results

command line option

```
--fps <fps>
                        moseq2-extract-convert-raw-to-avi --threads <threads>
                                                                                moseq2-extract-convert-raw-to-avi
                        command line option
                                                                                command line option
                        moseq2-extract-copy-slice
                                                                                moseq2-extract-copy-slice
                        command line option
                                                                                command line option
                        moseq2-extract-extract commandking-model-init <tracking_model_init>
                                                                                                  moseq2-extract-ex
                                                                                                  command line opt
                        line option
                                     moseq2-extrack-extrackodel-ll-clip <tracking_model_ll_clip>
                                                                                                      moseq2-extra
    --frame-dtype <frame_dtype>
                                                                                                      command line
                                     command line option
                                    moser2ckingtmextact-ll-threshold <tracking_model_ll_threshold>
                                                                                                           moseq
     --frame-trim <frame_trim>
                                                                                                           comma
                                    command line option
                               moseq2textractinggregate-results-threshold <tracking_model_mask_threshold>
                                                                                                              mos
  --input-dir <input_dir>
                                                                                                              com
                               command line option
                              moseq2-extract-generate-index
                                                                                                      moseq2-extra
                                                                                                      command line
                               command line option
                                                                                    moseq2-extract-extract
                                                            --use-plane-bground
     --max-height <max_height>
                                    moseq2-extract-extract
                                                                                    command line option
                                    command line option
                                                                                    moseq2-extract-find-roi
                                    moseq2-extract-extract
     --min-height <min height>
                                                                                    command line option
                                    command line option
moseq2-extract-ext
                                                                                                 command line optic
                                                 command line option
                                moseq2-extract-aggregate-resMerbose <verbose>
                                                                                moseq2-extract-convert-raw-to-avi
 --output-dir <output_dir>
                                                                                command line option
                                command line option
                                                                                moseq2-extract-extract command
                                moseq2-extract-download-flip-file
                                                                               line option
                                command line option
                                moseq2-extract-extract command
--write-movie <write_movie>
                                                                                         moseq2-extract-extract
                                                                                         command line option
                                line option
                                                                              moseq2-extract-generate-index
                                moseq2-extract-find-roi command
                                                                              command line option
                                line option
                                                                           moseq2-extract-convert-raw-to-avi
                                 moseq2-extract-convert-raw-to-avi
--output-file <output file>
                                                                           command line option
                                 command line option
                                                                           moseq2-extract-copy-slice
                                 moseq2-extract-copy-slice
                                                                           command line option
                                 command line option
                                                            -c
                                                                                  moseq2-extract-copy-slice
                                 moseq2-extract-generate-config
                                                                                  command line option
                                 command line option
                                                                                  moseq2-extract-extract
                                 moseq2-extract-generate-index
                                                                                  command line option
                                 command line option
                                                            -f
                                                                            moseq2-extract-aggregate-results
                              moseq2-extract-generate-index
                                                                            command line option
    --pca-file <pca_file>
                              command line option
                                                                            moseg2-extract-generate-index
atial-filter-size <spatial_filter_size>
                                              moseq2-extract-extract
                                                                            command line option
                                              command line option
                                                                            moseq2-extract-aggregate-results
                             moseq2-extract-generate-index
     --subpath <subpath>
                                                                            command line option
                             command line option
                                                                            moseq2-extract-generate-index
                                            moseq2-extract-extract
tail-filter-iters <tail_filter_iters>
                                                                            command line option
                                            command line option
tail-filter-shape <tail_filter_shape>
                                            moseq2-extract-extract
                                            command line option
                                           moseq2-extract-extract
-tail-filter-size <tail_filter_size>
                                           command line option
```

moseq2-extract-extract command line option

poral-filter-size <temporal_filter_size>

-0	moseq2-extract-aggregate-results command line option	convert_mkv_to_avi() (in module moseq2_extract.io.video)		
moseq2-extract-convert-raw-to-avi command line option		convert_pxs_to_mm() (in module moseq2_extract.util)		
		convert_raw_to_avi_function() (in module		
	moseq2-extract-copy-slice command line option	<pre>moseq2_extract.util) copy_h5_metadata_to_yaml_wrapper() (in module</pre>		
	moseq2-extract-generate-config command line option	moseq2_extract.helpers.wrappers)		
	moseq2-extract-generate-index command line option	moseq2_extract.helpers.data)		
-p	moseq2-extract-generate-index command line option	create_extract_h5() (in module moseq2_extract.helpers.data)		
-s	moseq2-extract-extract command line option	crop_and_rotate_frames() (in module moseq2_extract.extract.proc)		
-t	moseq2-extract-convert-raw-to-avi command line option	D		
	moseq2-extract-copy-slice	dict_to_h5() (in module moseq2_extract.util)		
	command line option	download_flip_command() (in module		
	moseq2-extract-extract command line option	moseq2_extract.gui)		
-v	moseq2-extract-convert-raw-to-avi	E		
	command line option	em_get_ll() (in module moseq2_extract.extract.track)		
Α		em_init() (in module moseq2_extract.extract.track)		
	results_command() (in module	em_iter() (in module moseq2_extract.extract.track)		
moseq2_extract.gu	i)	em_tracking() (in module moseq2_extract.extract.track)		
aggregate_extract_moseq2_extract.he	* * * * * * * * * * * * * * * * * * * *	escape_path() (in module moseq2_extract.util)		
apply_roi() (in module moseq2_extract.extract.proc)		extract_chunk() (in module moseq2_extract.extract.extract)		
		extract_command() (in module moseq2_extract.gui)		
B		extract_found_sessions() (in module moseq2_extract.gui)		
<pre>build_manifest() moseq2_extract.he</pre>	(in module lpers.data)	extract_wrapper() (in module		
build_path() (in mo	dule moseq2_extract.util)	moseq2_extract.helpers.wrappers)		
С		F		
	(in module moseq2_extract.util)	feature_hampel_filter() (in module moseq2_extract.extract.proc)		
check_progress() (in module moseq2_extract.gui) clean_dict() (in module moseq2_extract.util)		find_roi_command() (in module moseq2_extract.gui)		
		flip_file_wrapper() (in module		
,	module moseq2_extract.util)	moseq2_extract.helpers.wrappers)		
clean_frames() moseq2_extract.ex	(in module tract.proc)	G		
click_param_annot	() (in module moseq2_extract.util)	<pre>gen_batch_sequence() (in module moseq2_extract.util)</pre>		
command_with_commoseq2_extract.uti	• •	generate_config_command() (in module moseq2_extract.gui)		
compute_scalars() moseq2_extract.ex	(in module tract.proc)	<pre>generate_index_command()</pre>		
CONFIG_FILE moseq2-extrac option	ct-download-flip-file command line	generate_index_wrapper() (in module moseq2_extract.helpers.wrappers)		

```
get_bbox() (in module moseq2_extract.extract.proc)
                                                           М
get_bground_im()
                               (in
                                               module
                                                           make_gradient() (in module moseq2_extract.util)
moseq2_extract.extract.proc)
                                                           model smoother()
                                                                                                          module
                                                                                          (in
get_bground_im_file()
                                 (in
                                               module
                                                           moseq2_extract.extract.proc)
moseq2_extract.extract.proc)
                                                           module
get_bucket_center() (in module moseq2_extract.util)
                                                               moseq2_extract.extract.extract
get_flips() (in module moseq2_extract.extract.proc)
                                                               moseg2 extract.extract.proc
get_found_sessions() (in module moseq2_extract.gui)
                                                               moseq2_extract.extract.roi
get_frame_features()
                                               module
                                (in
                                                               moseq2 extract.extract.track
moseg2 extract.extract.proc)
                                                               moseq2_extract.gui
                              (in
                                               module
get_largest_cc()
                                                               moseq2_extract.helpers.data
moseq2_extract.extract.proc)
                                                               moseg2 extract.helpers.extract
get_movie_info() (in module moseq2_extract.io.video)
                                                               moseq2_extract.helpers.wrappers
get_raw_info() (in module moseq2_extract.io.video)
                                                               moseq2 extract.io.image
get_roi() (in module moseq2_extract.extract.proc)
                                                               moseq2_extract.io.video
get_roi_wrapper()
                               (in
                                               module
moseq2_extract.helpers.wrappers)
                                                               moseq2_extract.util
get selected sessions()
                                  (in
                                               module
                                                           moseq2-extract-aggregate-results command line
moseq2_extract.helpers.data)
                                                           option
                                                               --format <format>
get_video_info() (in module moseq2_extract.io.video)
                                                               --input-dir <input_dir>
graduate_dilated_wall_area()
                                    (in
                                               module
moseq2_extract.util)
                                                               --output-dir <output dir>
                                                               -f
H
                                                               -i
h5 to dict() (in module moseg2 extract.util)
handle_extract_metadata()
                                   (in
                                               module
                                                           moseq2-extract-convert-raw-to-avi command line
moseq2_extract.helpers.data)
                                                           option
                                                               --chunk-size <chunk_size>
                                                               --delete
                                 (in
                                               module
im moment features()
                                                               --fps <fps>
moseq2 extract.extract.proc)
INPUT_FILE
                                                               --output-file <output_file>
    moseg2-extract-convert-raw-to-avi command line
                                                               --threads <threads>
    option
                                                               --verbose <verbose>
    moseq2-extract-copy-slice command line option
                                                               -b
    moseg2-extract-extract command line option
                                                               -0
    moseq2-extract-find-roi command line option
                                                               -t
L
                                                               -v
                                                               INPUT FILE
load_h5s() (in module moseq2_extract.helpers.data)
                                                           moseq2-extract-copy-slice command line option
load metadata() (in module moseg2 extract.util)
                                                               --chunk-size <chunk_size>
load movie data()
                                               module
                               (in
                                                               --copy-slice <copy_slice>
moseq2 extract.io.video)
                                                               --delete
load_textdata() (in module moseq2_extract.util)
                                                               --fps <fps>
load_timestamps() (in module moseq2_extract.util)
                                                               --output-file <output_file>
                                                               --threads <threads>
```

-b --frame-trim <frame trim> --max-height <max_height> -C --min-height <min_height> -0 --model-smoothing-clips <model_smoothing_clips> -t INPUT_FILE --output-dir <output_dir> moseg2-extract-download-flip-file command line --spatial-filter-size <spatial_filter_size> option --tail-filter-iters <tail filter iters> --output-dir <output_dir> --tail-filter-shape <tail_filter_shape> CONFIG FILE --tail-filter-size <tail_filter_size> moseq2-extract-extract command line option --temporal_filter_size <temporal_filter_size> --angle-hampel-sig <angle_hampel_sig> --tracking-model-init <tracking_model_init> --angle-hampel-span <angle_hampel_span> --tracking-model_ll_clip <tracking_model_ll_clip> --bg-roi-depth-range <bg_roi_depth_range> --tracking-model-II-threshold --bg-roi-dilate <bg_roi_dilate> <tracking model II threshold> --bg-roi-fill-holes <bg_roi_fill_holes> --tracking-model-mask-threshold --bg-roi-gradient-filter <bg_roi_gradient_filter> <tracking_model_mask_threshold> --bg-roi-gradient-kernel <bg_roi_gradient_kernel> --tracking-model-segment <tracking_model_segment> --bg-roi-gradient-threshold <bg_roi_gradient_threshold> --use-plane-bground --bg-roi-index <bg_roi_index> --use-tracking-model <use_tracking_model> --bg-roi-shape <bg_roi_shape> --verbose <verbose> --bg-roi-weights <bg_roi_weights> --write-movie <write movie> --bg-sort-roi-by-position <bg_sort_roi_by_position> -C --bg-sort-roi-by-position-max-rois -s <bg_sort_roi_by_position_max_rois> -t --cable-filter-iters <cable_filter_iters> INPUT_FILE --cable-filter-shape <cable_filter_shape> moseq2-extract-find-roi command line option --cable-filter-size <cable filter size> --bg-roi-depth-range <bg_roi_depth_range> --centroid-hampel-sig <centroid_hampel_sig> --bg-roi-dilate <bg roi dilate> --centroid-hampel-span <centroid_hampel_span> --bg-roi-fill-holes <bg_roi_fill_holes> --chunk-overlap <chunk overlap> --bg-roi-gradient-filter <bg roi gradient filter> --chunk-size <chunk_size> --bg-roi-gradient-kernel <bg roi gradient kernel> --compress <compress> --bg-roi-gradient-threshold <bg_roi_gradient_threshold> --compress-chunk-size <compress_chunk_size> --bg-roi-index <bg_roi_index> --compress-threads <compress_threads> --bg-roi-shape <bg roi shape> --config-file <config_file> --bg-roi-weights <bg_roi_weights> --crop-size <crop_size> --bg-sort-roi-by-position <bg_sort_roi_by_position> --detected-true-depth <detected_true_depth> --bg-sort-roi-by-position-max-rois --dilate-iterations <dilate_iterations> <bg_sort_roi_by_position_max_rois> --flip-classifier <flip_classifier> --config-file <config_file> --flip-classifier-smoothing --dilate-iterations <dilate iterations> <flip_classifier_smoothing> --output-dir <output dir> --fps <fps> --use-plane-bground --frame-dtype <frame_dtype>

```
moseq2_extract.helpers.extract)
moseq2-extract-generate-config command line
option
                                                           R
    --output-file <output file>
                                                          read_frames() (in module moseq2_extract.io.video)
    -0
                                                                                          (in
                                                          read frames raw()
                                                                                                         module
moseq2-extract-generate-index command line
                                                          moseg2 extract.io.video)
option
    --all-uuids <all_uuids>
                                                          read_image() (in module moseq2_extract.io.image)
    --filter <filter>
                                                          read_yaml() (in module moseq2_extract.util)
    --input-dir <input dir>
                                                          recursive find h5s() (in module moseg2 extract.util)
    --output-file <output_file>
                                                          recursive_find_unextracted_dirs()
                                                                                                 (in
                                                                                                         module
                                                          moseq2_extract.util)
    --pca-file <pca_file>
                                                                                            (in
                                                                                                         module
                                                          restore_progress_vars()
    --subpath <subpath>
                                                          moseq2_extract.gui)
    -a
                                                          run local extract()
                                                                                          (in
                                                                                                         module
    -f
                                                          moseq2_extract.helpers.extract)
    -i
                                                          run slurm extract()
                                                                                                         module
                                                                                          (in
                                                          moseq2 extract.helpers.extract)
    -0
    -p
                                                           S
moseq2_extract.extract.extract
                                                                                              (in
                                                                                                         module
                                                          sample_extract_command()
    module
                                                          moseq2_extract.gui)
moseq2_extract.extract.proc
                                                          scalar attributes() (in module moseq2 extract.util)
    module
                                                          select_strel() (in module moseq2_extract.util)
moseq2_extract.extract.roi
                                                          strided_app() (in module moseq2_extract.util)
    module
moseq2_extract.extract.track
                                                           T
    module
moseq2 extract.gui
                                                          time_str_for_filename() (in module moseq2_extract.util)
    module
                                                           U
moseq2_extract.helpers.data
    module
                                                          update_progress() (in module moseq2_extract.gui)
moseq2_extract.helpers.extract
    module
moseq2_extract.helpers.wrappers
                                                          view_extraction() (in module moseq2_extract.gui)
    module
moseq2_extract.io.image
                                                           W
    module
                                                          write_frames() (in module moseq2_extract.io.video)
moseq2 extract.io.video
                                                                                                         module
                                                          write frames preview()
    module
                                                          moseq2_extract.io.video)
moseq2 extract.util
                                                          write_image() (in module moseq2_extract.io.image)
    module
mouse threshold filter()
                                 (in
                                              module
moseg2 extract.util)
```

INPUT FILE

plane_fit3() (in module moseq2_extract.extract.roi)
plane_ransac() (in module moseq2_extract.extract.roi)

(in

process_extract_batches()

module

Python Module Index

m

moseq2_extract
moseq2_extract.extract.extract.extract
moseq2_extract.extract.proc
moseq2_extract.extract.roi
moseq2_extract.extract.track
moseq2_extract.gui
moseq2_extract.helpers.data
moseq2_extract.helpers.extract
moseq2_extract.helpers.wrappers
moseq2_extract.io.image
moseq2_extract.io.video

moseq2_extract.util